

How Do I Choose the Right Microchip Transplant?

Features	Datamars	Others
1. Reputation – <i>Is it an internationally recognized product and certified by ICAR and ISO?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
2. ISO Standard - <i>Does the product comply with ISO 11784 & 11785 with worldwide compatibility?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
3. Coding – <i>Does the transponder provide manufacture and country coding?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
4. Sterility Date – <i>Is the product single packed with sterilized coloured indication and sterility expiry date?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
5. Disposal Packing – <i>Does it come with one-way disposal syringe for one application to prevent spreading of diseases / infections?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
6. Reading Distance – <i>Is the reading distance sufficient to ensure fast and reliable detection by the scanner?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
7. Needle Length – <i>Is the needle length of appropriate measurement (about 25mm max in length) in order not to penetrate to deep and not to injure the animal?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
8. Ergonomic Design – <i>Is the needle ergonomically designed for better finger grip control and easy handling?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
9. Safety Feature – <i>Does the syringe come with a plunger stopper (click feature) to ensure proper placement of microchip once implanted?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
10. ID Label – <i>Does the transponder pack contain 6 adhesive labels with ID numbers?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
	TOTAL	10

Peter Huber
Director
Efre Far East Pte Ltd