

## How Do I Choose the Right Microchip Transplant?

Features	Datamars	Others
1. <b>Reputation</b> – <i>Is it an internationally recognized product and certified by ICAR and ISO?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
2. <b>ISO Standard</b> - <i>Does the product comply with ISO 11784 &amp; 11785 with worldwide compatibility?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
3. <b>Coding</b> – <i>Does the transponder provide manufacture and country coding?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
4. <b>Sterility Date</b> – <i>Is the product single packed with sterilized coloured indication and sterility expiry date?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
5. <b>Disposal Packing</b> – <i>Does it come with one-way disposal syringe for one application to prevent spreading of diseases / infections?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
6. <b>Reading Distance</b> – <i>Is the reading distance sufficient to ensure fast and reliable detection by the scanner?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
7. <b>Needle Length</b> – <i>Is the needle length of appropriate measurement (about 25mm max in length) in order not to penetrate to deep and not to injure the animal?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
8. <b>Ergonomic Design</b> – <i>Is the needle ergonomically designed for better finger grip control and easy handling?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
9. <b>Safety Feature</b> – <i>Does the syringe come with a plunger stopper (click feature) to ensure proper placement of microchip once implanted?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
10. <b>ID Label</b> – <i>Does the transponder pack contain 6 adhesive labels with ID numbers?</i>	<input checked="" type="checkbox"/>	<input type="checkbox"/>
<b>TOTAL</b>	<b>10</b>	

Peter Huber  
Director  
Efre Far East Pte Ltd